

CHOOSING THE RIGHT FINISH FOR YOUR APPLICATION

When our engineers receive enquiries for bespoke fabrication, one of the first questions they will ask is “what finish do you require?”.

This is because the finish can greatly effect the cost, quality and in some cases delivery times. Not every enquiry requires a pharmaceutical finish with full certification and testing!

Axiom’s range of finishes include:


Orbital: Welds left as laid, no post weld treatment.

Manual: Welds left as laid, pickled and brushed.


0.8Ra: Welds removed internally and externally, polished 0.8Ra.

0.4Ra: Welds removed internally and externally, polished to 0.4Ra (top end pharmaceutical, high value biotechnology components).

Mirror polish: Welds removed internally and externally, Mirror polish (top end pharmaceutical, high value biotechnology components).



ORBITAL



MANUAL



0.8Ra



0.4Ra



MIRROR POLISH



www.axiumprocess.com

Henly Industrial Estate, Henly, Swansea SA4 0XP UK

Phone: +44 (0)1792 883882 **Fax:** +44 (0)1792 886049

Email: info@axiumprocess.com

Axiom Process Ltd; Company registered in England & Wales No.: 4432600; VAT registration: 793 7784 59



FS 83081